

# *Celloger<sup>®</sup> Stack*

*Automated multi-layer  
vessel monitoring system*



*“ Find out the right time to harvest your cells  
without opening the incubator ”*





## Key features

- » Real-time live cell monitoring & multi-layer observation (up to 10 layers)
- » Multi-position time-lapse imaging & recording
- » Checks cell quality and detects confluency rate while culturing
- » Improved reproducibility for mass cell culture and production
- » No need to take the cells out of the incubator thus eradicating contamination worries

# CURIOSIS

## Cell images taken by the Celloger® Stack


Hela cell incubated in a 10-layer chamber for 7 days with 1 hour interval time



## Specification

Dimensions	350 x 330 x 450 mm
Weight	15 kg
Objective Lens	2X
Camera	5MP CMOS
Stage	Motorized XYZ stage
Operating environment	10~40°C, 20~95% humidity
File export format	TIFF, AVI
Resolution	2592 x 1944 pixel
Imaging positions	Multiple
Imaging mode	Bright-field
Vessel types	Multi-layer chamber up to 10 layers
Key functions	Autofocusing, manual focusing, time-lapse imaging, movie maker, confluency, growth curve

## Confluency Graph



Contact us for more information about the system

12086 Miramar Pkwy  
Miramar, Florida, 33015, USA  
www.prolabcorp.com  
+1 305 6000 444

