

Conical Tube (#50715/50750)

1. Description

SPL Life Sciences provides a new type of Conical Tube including highly resilient silicone septum in the center of the cap. Septum Tube is effective in completely blocking external contaminants because it is not necessary to open the cap to handle the sample. Silicone septum with PTFE coating is highly resistant to various solvents and fully recovers even repeated use of syringes.

2. Features

- ❖ Tight sealing screw cap
- ❖ Silicone septum with PTFE coating
- ❖ External graduations with marking area
- ❖ Packed in zipper bags
- ❖ Provided in racks

3. Functional Test

Non-Cytotoxic : Passed
 Non-Pyrogenic : Passed
 DNase/RNase Free : Passed
 Human DNA Free : Passed

4. Product Range

Cat.No.	Material (Tube/Cap)	Septum Material	External Dimensions d x h (mm)	Total Vol. (ml)	RCF Rating	Sterile	Packaging
50715	PP / HDPE	Silicone / PTFE	17.00 x 120.00	15.00	13,000 x g	+	25 in rack
50750	PP / HDPE	Silicone / PTFE	30.00 x 115.00	50.00	14,000 x g	+	25 in rack

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs



SPL-TDS-COTU


7. Specifications & Technical Drawing

(1) 15 ml Septum Tube (# 50715)

- Specifications

	Dimensions		Material (Resin)	Packaging
Cap	Diameter:	21.34 mm	High-density polyethylene (HDPE)	• Pieces/Case: 25
	Height:	13.30 mm		
Body	Diameter:	16.90 mm	Polypropylene(PP)	
	Height:	119.20 mm		
	Volume:	15.00 ml		


- Technical Drawing



(2) 50 ml Septum Tube (# 50750)

- Specifications

	Dimensions		Material (Resin)	Packaging
Cap	Diameter:	35.02 mm	High-density polyethylene (HDPE)	• Pieces/Case: 25
	Height:	18.50 mm		
Body	Diameter:	29.30 mm	Polypropylene(PP)	
	Height:	113.90 mm		
	Volume:	50.00 ml		



12086 Miramar Pkwy
 Miramar, Florida, 33025, USA
 +1 (305) 600-0444
 info@prolabcorp.com
 www.prolabcorp.com