



Sample Tube (#50221)

1. Description

SPL Sample Tubes are designed for storing various clinical samples. SPL Life Sciences provides durable, easy-to-use bottles to meet various demands of customers.

2. Features

- ❖ Tight seal with screw cap
- ❖ External graduations
- ❖ Self-standing type conical bottom
- ❖ Autoclavable

3. Product Range

Cat. No.	Material (Body/Cap)	External Dimension (mm)	Total Volume (ml)	Sterile	Packaging
50221	PP / PE	25.00 x 73.00	20.00	-	25 / 500

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



7. Specifications & Technical Drawing

- Specifications

Dimensions	Material (Resin)	Packaging
Diameter: 25.20 mm Height: 75.00 mm	<ul style="list-style-type: none"> • Polypropylene (PP, body) • High Density Polyethylene (HDPE, Cap) 	<ul style="list-style-type: none"> • Pieces/Sleeve: 25 • Pieces/Case: 500

- Technical Drawing

