

Wide-Mouth Bottle (PP)

(#58005/58020/58030/58031/58060/58125/58250/58500/58000)

1. Description

SPL Wide-Mouth Bottles are durable bottles for storing and handling liquid/solid forms for samples and chemicals. Various sizes are available.

2. Features

- ❖ Volumes range from 5 to 1,000 ml
- ❖ 250 ml bottles are applicable for centrifugation (RCF 4,000 g)
- ❖ Colored cap available for 5 ml bottles : Blue, Green, Red, Yellow (Cat. No. 58005)
- ❖ Autoclavable

3. Product Range

Cat. No.	Material (Bottle / Cap)	Total Volume (ml)	Packaging
58005	PP / PP	5.00	1000
58020	PP / PP	20.00	1000
58030	PP / PP	30.00	1000
58031	PP / PP	30.00	12 / 72
58060	PP / PP	60.00	1000
58125	PP / PP	125.00	500
58250	PP / PP	250.00	72
58500	PP / PP	500.00	48
58000	PP / PP	1,000.00	6 / 24

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs




SPL-TDS-WMVPP

6. Specifications & Technical Drawing


- Specifications

	Dimensions			Material (Resin)	Packaging	
	58005	58020	58030/58031		58005/58020/58030	58031
Diameter	18.00 mm	32.00 mm	34.20 mm	•Body: PP •Cap: PP	• Pieces/Case: 1000	• Pieces/Sleeve: 12 • Pieces/Case: 72
High	42.00 mm	53.50 mm	61.80 mm			


- Technical Drawing



Cat. No. 58005



Cat. No. 58020




Cat. No. 58030/58031


- Specifications

	Dimensions				Material (Resin)	Packaging			
	58060	58125	58250	58500		58060	58125	58250	58500
Diameter	38.50 Ø	50.20 Ø	61.00 Ø	72.00 Ø	•Body: PP •Cap: PP	Pieces/Case: 1000	Pieces/Case: 500	Pieces/Case: 72	Pieces/ Case: 48
High	84.80mm	97.30mm	130.00mm	166.30 mm					


- Technical Drawing




Cat. No. 58060



Cat. No. 58125



Cat. No. 58250




Cat. No. 58500

- Specifications

	Dimensions	Material (Resin)	Packaging
	58000		58000
Diameter	91.00 mm	•Body: PP •Cap: PP	• Pieces/Case: 24
High	196.00 mm		

- Technical Drawing



Cat. No. 58000